MULTI SCIENCES

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Mouse Regulatory T Cell Staining Kit

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This package insert must be read entirely before using this product. For proper performance, use the insert provided with each individual product received.

Catalog Number

KTR201 - 25

KTR201 - 100

Optimization for mouse regulatory T cell (Treg) staining in anticoagulated blood and splenocytes.

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For research use only. Not for use in diagnostic procedures.

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INTRODUCTION

The regulatory T cells (Tregs), formerly known as suppressor T cells, are a subpopulation of T cells that modulate the immune system, maintain tolerance to self-antigens, and prevent autoimmune disease. Tregs are immunosuppressive and generally suppress or downregulate induction and proliferation of effector T cells. Regulatory T cells come in many forms with the most well-understood being those that express CD4, CD25, and FOXP3 (CD4⁺CD25⁺ regulatory T cells).

The immune system must be able to discriminate between self and non-self. When self/non-self discrimination fails, the immune system destroys cells and tissues of the body and as a result causes autoimmune diseases. Regulatory T cells actively suppress activation of the immune system and prevent pathological self-reactivity, i.e. autoimmune disease. The critical role regulatory T cells play within the immune system is evidenced by the severe autoimmune syndrome that results from a genetic deficiency in regulatory T cells (IPEX syndrome).

Mouse models have suggested that modulation of Tregs can treat autoimmune disease and cancer and can facilitate organ transplantation. Their implications for cancer are complicated. Tregs tend to be upregulated in individuals with cancer, and they seem to be recruited to the site of many tumors. Studies in both humans and animal models have implicated that high numbers of Tregs in the tumor microenvironment is indicative of a poor prognosis, and Tregs are thought to suppress tumor immunity, thus hindering the body's innate ability to control the growth of cancerous cells. Recent immunotherapy research is studying how regulation of T cells could possibly be utilized in the treatment of cancer.

LIMITATIONS OF THE PRODUCT

The product is intended for flowcytometry applications. FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The product should not be used beyond the expiration date on the label.

Do not mix reagents with those from other lots or sources.

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

MATERIALS PROVIDED AND STORAGE

Components	Product Code	-25	-100
Anti-Mouse CD4, FITC (Clone: GK1.5)	AM00401	150 µl	600 µl
Anti-Mouse CD25, APC (Clone: PC61.5)	AM02505	150 µl	600 µl
Anti-Mouse FoxP3, PE (Clone: 3G3)	AM0F04	150 µl	600 µl
Anti-Mouse CD16/CD32, Purified (Clone: 2.4G2)	AM016	60 µl	240 µl
Mouse IgG1 Isotype Control, PE	CMG104	25 µl	100 µl
Fixation/Permeabilization Concentrate (4×)	ICC01	8 ml	30 ml
Fixation/Permeabilization Diluent	ICD01	25 ml	100 ml
Permeabilization Buffer (10×)	ICB01	25 ml	100 ml
FCM Lysing Solution (Fixative Free) (10×)	LYS	8 ml	30 ml
Flow Cytometry Staining Buffer (1×)	S1001	125 ml	125 ml × 3

Note: All reagents can be stored at 2 - 8°C, and would be stable for at least 1 year when stored at recommended condition.

OTHER SUPPLIES REQUIRED

1. **Distilled water.**
2. **Rat IgG2b Isotype Control, FITC** (Cat No. CRG2b01, MultiSciences) (Isotype control for anti-mouse CD4, FITC, Maybe needed)
3. **Rat IgG1 Isotype Control, APC** (Isotype control for anti-mouse CD25, APC, Maybe needed)
4. **12 × 75 mm round bottom test tubes.**
5. **Vortexer.**
6. **Swing-out horizon centrifuge** (with rotor for 15 ml tubes).

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PROTOCOL

Sample collection and treatment

Collect appropriate amount of anticoagulant blood, store at room temperature or 2 - 8°C, and detect on the same day.

For spleen, collect and prepare as single cell suspension by using appropriate method, and remove tissue mass, detect on the same day.

Reagent preparation


Prepare 1× FCM Lysing Solution (Fixation Free) by diluting the 10× concentrate with distilled water prior to use. You will need 2 ml of 1× FCM Lysing Solution for each sample, if staining in tubes. Store excess at 2 - 8°C.

Prepare fresh Fixation/Permeabilization working solution by diluting the Fixation/Permeabilization Concentrate (1 part) with Fixation/Permeabilization Diluent (3 parts). You will need 1 ml of the Fixation/Permeabilization working solution for each sample, if staining in tubes. Do not store this buffer more than 1 day.

Prepare a 1× working solution of Permeabilization Buffer by diluting the 10× concentrate with distilled water prior to use. You will need 8.5 ml of Permeabilization Buffer for each sample, if staining in tubes. Store excess at 2 - 8°C for up to 1 week.

Sample preparation

1. Pipet 100 µl anticoagulant blood or $1 - 10 \times 10^6$ splenocytes into test tubes, add 5 µl Anti-Mouse CD4, FITC and 5 µl Anti-Mouse CD25, APC to each tube. Vortex to mix well, and incubate at room temperature for 15 minutes, protect from light.
2. Add 2 ml of 1× FCM Lysing Solution work solution to each tube and pulse vortex, and incubate at room temperature for 15 minutes, protect from light. Centrifuge at $300 - 400 \times g$ for 5 minutes at room temperature, then discard the supernatant. Add 2 ml 1× Flow Cytometry Staining Buffer, and pulse vortex. Centrifuge at $300 - 400 \times g$ for 5 minutes at room temperature, then discard the supernatant.
3. Add 1 ml of Fixation/Permeabilization working solution to each tube and pulse vortex.

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
4. Incubate for 30 - 60 minutes at room temperature. Protect samples from light.
5. Without washing, add 2 ml of 1× Permeabilization Buffer to each tube. Centrifuge at 300 - 400 × g for 5 minutes at room temperature, then discard the supernatant.
6. (Optional) Repeat step 5.
7. Resuspend pellet in 100 µl of 1× Permeabilization Buffer. This is typically the residual volume after decanting.
8. (Optional) Block with Anti-Mouse CD16/CD32, Purified by adding 2 µl directly to the cells. Incubate for 15 minutes at room temperature, protect from light.
9. Without washing, add 5 µl Anti-Mouse Foxp3, PE or 5 µl Mouse IgG1 Isotype Control, PE to cells. Vortex to mix well, and incubate for at least 30 minutes at room temperature and protect samples from light.
10. Add 2 ml of 1× Permeabilization Buffer to each tube. Centrifuge at 300 - 400 × g for 5 minutes at room temperature, then discard the supernatant.
11. (Optional) Repeat step 10.
12. Resuspend stained cells in 500 µl 1× Flow Cytometry Staining Buffer and acquire samples on a flow cytometer.


Detection by flowcytometry

Correct gating to determine the frequencies of Tregs in samples.

Note: Acquire at least 20,000 to 30,000 CD4⁺ cells. In order to make statistically significant frequency measurements, sufficiently large sample sizes should be acquired during flow cytometric analysis.

Excitation wavelength and emission wavelength of fluorophore

Fluorophore	Ex (nm)	Em (nm)
FITC	495	519
R-Phycoerythrin (PE)	480;565	578
APC	650	660

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

Examples of results/结果示例

Flow cytometric analysis of Mouse Regulatory T Cell Staining Kit. The staining pattern of FSC/CD4 (Left) CD25/Mouse IgG1 Isotype Control, PE (Middle) and CD25/FoxP3 (Right) on splenocytes of normal ICR mouse. Adherent cells have been excluded. Flow cytometry was performed on a Thermo Fisher Attune NxT.

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