



5× Binding Buffer

I. 产品信息

目 录 号: AP101-30-BB, AP101-60-BB, AP101-100-BB 规 格: 5ml, 10ml, 15ml
保 存: 2 - 8℃ 效 期: 一年

II. 背景简介

Annexin V(或 Annexin A5)为胞内蛋白膜联蛋白家族成员, 以钙依赖的方式与磷脂酰丝氨酸(PS)结合。PS 存在于正常细胞浆膜的内层, 但在凋亡早期, 膜不对称性丧失, PS 易位至细胞表面。

本试剂可提供荧光标记的 Annexin V 与 PS 结合的环境, 可与其他单独组分如 Annexin V-FITC、PI 等搭配使用, 也可购买完整的试剂盒。

III. 使用方法

使用前, 用双蒸水将 5 × Binding Buffer 稀释为 1 × 工作液。

IV. 注意事项

1. 请在使用本产品前仔细阅读说明书。本产品仅用于科研, 不可用于诊断。
2. 为了您的安全和健康, 请穿戴实验防护服、手套、口罩等必要的防护装备。
3. 更多凋亡相关产品敬请关注联科生物网站或来电咨询。

V. 部分相关产品

货号	产品名称	规格
AP101-30	Annexin V-FITC/PI apoptosis kit	30T
AP101-30-AVF	Annexin V-FITC	30T
AP101-30-PI	PI solution (For apoptosis detection with Annexin V)	30T
AP104-30	Annexin V-PE/7-AAD apoptosis kit	30T
AP104-30-AAD	7-AAD	30T
AP104-30-AVP	Annexin V-PE	30T
AP105-30-AVA	Annexin V-APC	30T
AP105-30	Annexin V-APC/7-AAD apoptosis kit	30T
AP107-30	Annexin V-APC/PI apoptosis kit	30T
APCC101-30	Annexin V-FITC/PI apoptosis and cell cycle kit	30T
APCC104-30	Annexin V-PE/7-AAD apoptosis and cell cycle kit	30T
APCC105-30	Annexin V-APC/7-AAD apoptosis and cell cycle kit	30T
C050	CCCP, 50mM	125 µl
CCS013	Apo-peak DNA analysis kit	50T
J200	JC-1, 2mM	125 µl
MJ101	Mitochondria Staining Kit (JC-1)	125T
00-4555-56	Accutase - Enzyme Cell Detachment Medium	100 ml

Product For Research Use Only rev: AP101-BB-01

MultiSciences Biotech Co., Ltd

Technical Service: tech@liankebio.com

To place an order: info@liankebio.com